

## CATHETERS FOR HSG/SIS

### FLEXIBLE HS CATHETER

- THIS PRODUCT DOES NOT CONTAIN LATEX
- Sterile. Disposable. Ready to use.
- Placement sheath decreases the need for a tenaculum
- Patient ambulatory once instrument properly inserted
- Allows for clear view of external cervical os
- Stylet available for cases when additional stiffness is needed
- For use with saline or water based contrast media
- 5fg and 7fg sizes available in a box of 10

### SHAPEABLE HS CATHETER

- THIS PRODUCT DOES NOT CONTAIN LATEX
- Sterile. Disposable. Ready to use.
- Integrated stylet is malleable and can be re-shaped as necessary
- Angle can be set to align with cervical canal
- Stiffness beneficial for restricted or stenotic cervical canal
- For use with saline or water-based contrast media
- 5fg and 7fg sizes available in a box of 10


*New "Twin Tap" HS Catheter now available*

## ExEm Foam-Kit

### Tubal patency with clear results

ExEm foam kits provide an inexpensive, simple, one kit solution for outpatient ultrasound testing of tubal patency and uterine anatomy in patients with fertility or gynaecological problems. Achieve better results in HyCosy procedures - the imaging is superior and longer lasting than any other sonohysterography products.

By mixing ExEm® gel and purified water a gel foam is created. After infusing the foam through the uterine cavity into the fallopian tubes, practical ultrasound images can be obtained to check the patency of the fallopian tubes in infertility patients. Subsequently high quality ultrasound images can be obtained. The ExEm® Foam-kit contains a 10 ml syringe with ExEm® gel, a 10 ml syringe with purified water, an empty 20 ml syringe and a coupling device.



# Ring Pessaries

The Wallace ring pessaries are available in a wide range of sizes to suit all patients. The flexible material provides comfort to the patient during fitting and duration of use.

Designed to ease uterine prolapse and other uterine displacements.

Simple to fit: Ring softens when immersed in warm water to aid fitting in a figure of eight configuration.

Easy to maintain: Remove and replace with a new ring pessary after three months.

## KEY PRODUCT BENEFITS

- Range of sizes to suit all patients
- Flexible PVC with large cross section
- Sixteen sizes from 50mm to 110mm (outside diameter)
- Individually packed in easy to open pouches

## PIPELLE DE CORNIER®

- Systematic screening for endometrial cancer and hyperplasia
- Detection of luteal phase insufficiency
- Monitoring of hormonal treatment
- Menometrorrhagia with or without HRT
- Late menopause
- Endometrial hypertrophy detected by ultrasonography
- Polyps
- Tamoxifen treatment
- Bacteriological culture

## Pipelle de Cornier® Mark II



The Pipelle de Cornier® Mark II is a new device based on the original Pipelle de Cornier®. The endometrial cells (cytology specimen) are collected on the rough surface on 2 cm of its tip, and the endometrial tissue (histology specimen) is collected by the plunger's aspiration.



## GYTECH DISPOSABLE UTERINE SOUND

- Sterile. Disposable, Ready to use
- Length 24.8cm
- Shapeable Tip
- OD 3.5mm
- This product does not contain latex
- Box of 25

## GYTECH FLEXIBLE CERVICAL DILATOR (OS FINDER)

- Sterile. Disposable, Ready to use
- Length 21.5cm
- OD 2-4mm tapered
- This product does not contain latex
- Phthalate Free
- Box of 25 or 50